

Neop | 600 x 400 distribution trolleys

Daily distribution of medicines with pill dispensers and nominative trays

- Robustness: excellent resistance, one piece structure made of aluminium
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ward name identification
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours



Robust



Hygienic



Ergonomic



Customizable

- Single column trolley
- 3 height: 80, 100 or 120 cm
- Roller shutter door or door wing
- Code lock (with spare key) or badge lock
- 4 available inside panel models (slides) according to internal fittings



Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue

Standard configurations - PREMIUM

Structure WITH internal fittings and peripheral accessories.

Rails position:

- Top right
- Bottom right
- Bottom left

Trolleys height 100 cm

DIN – 18 patients

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 18 bins H70 l.1/3 – white color
- 7 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



| Reference | Description |
|-------------------|-----------------------------------------------------------------------------|
| 8YM64H100DIN18_RM | NEOP "Premium" medicine trolley 600x400 H100 DIN 18 patients – magenta pink |
| 8YM64H100DIN18_BC | NEOP "Premium" medicine trolley 600x400 H100 DIN 18 patients – sky blue |
| 8YM64H100DIN18_VP | NEOP "Premium" medicine trolley 600x400 H100 DIN 18 patients – apple green |

Label covers to order if necessary

DIN – 24 patients

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 24 bins H70 L.1/4 – white color
- 7 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H75
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



| Reference | Description |
|-------------------|-----------------------------------------------------------------------------|
| 8YM64H100DIN24_RM | NEOP "Premium" medicine trolley 600x400 H100 DIN 24 patients – magenta pink |
| 8YM64H100DIN24_BC | NEOP "Premium" medicine trolley 600x400 H100 DIN 24 patients – sky blue |
| 8YM64H100DIN24_VP | NEOP "Premium" medicine trolley 600x400 H100 DIN 24 patients – apple green |

Label covers to order if necessary

DIN – 28 patients

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 28 bins H70 L.1/4 – white color
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket




| Reference | Description |
|-------------------|----------------------------------------------------------------------------|
| 8YM64H100DIN28_RM | NEOP "Premium" medicine trolley 600x400 H100 DIN 28 lits – magenta pink |
| 8YM64H100DIN28_BC | NEOP "Premium" medicine trolley 600x400 H100 DIN 28 lits – sky blue |
| 8YM64H100DIN28_VP | NEOP "Premium" medicine trolley 600x400 H100 DIN 28 patients – apple green |

Label covers to order if necessary

PDA PILI – 70 patients

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 batches of 10 PILI dispenser supports
- 7 trays 600 x 400 for PILI delivered without dispenser
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket




| Reference | Description |
|-------------------|----------------------------------------------------------------------------------|
| 8YM64H100PIL70_RM | NEOP "Premium" medicine trolley 600x400 H100 PDA Pili 70 patients - magenta pink |
| 8YM64H100PIL70_BC | NEOP "Premium" medicine trolley 600x400 H100 PDA Pili 70 patients - sky blue |
| 8YM64H100PIL70_VP | NEOP "Premium" medicine trolley 600x400 H100 PDA Pili 70 patients - apple green |

Label covers to order if necessary

PDA MODULO CLASSIC – 70 patients

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 7 trays 600 x 400 with 10 Modulo Classic dispensers
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



| Reference | Description |
|------------------|--------------------------------------------------------------------------------------------|
| 8YM64H100MC70_RM | NEOP "Premium" medicine trolley 600x400 H100 PDA Modulo Classic 70 patients - magenta pink |
| 8YM64H100MC70_BC | NEOP "Premium" medicine trolley 600x400 H100 PDA Modulo Classic 70 patients - sky blue |
| 8YM64H100MC70_VP | NEOP "Premium" medicine trolley 600x400 H100 PDA Modulo Classic 70 patients - apple green |

Label covers to order if necessary

PDA ROLL DISPENSER – 24 patients

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 24 bins H110 L.1/6 + roll dispenser + divider
- 5 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H75
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket




| Reference | Description |
|-------------------|-------------------------------------------------------------------------------------------------|
| 8YM64H100DEV24_RM | NEOP "Premium" medicine trolley 600x400 H100 PDA roll dispenser H110 24 patients – magenta pink |
| 8YM64H100DEV24_BC | NEOP "Premium" medicine trolley 600x400 H100 PDA roll dispenser H110 24 patients – sky blue |
| 8YM64H100DEV24_VP | NEOP "Premium" medicine trolley 600x400 H100 PDA roll dispenser H110 24 patients – apple green |

Label covers to order if necessary

PDA ROLL DISPENSER – 30 patients

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 30 bins H110 L.1/6 + roll dispenser + divider
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



| Reference | Description |
|-------------------|-------------------------------------------------------------------------------------------------|
| 8YM64H100DEV30_RM | NEOP "Premium" medicine trolley 600x400 H100 PDA roll dispenser H110 30 patients - magenta pink |
| 8YM64H100DEV30_BC | NEOP "Premium" medicine trolley 600x400 H100 PDA roll dispenser H1100 30 patients - sky blue |
| 8YM64H100DEV30_VP | NEOP "Premium" medicine trolley 600x400 H100 PDA roll dispenser H110 30 patients - apple green |

Label covers to order if necessary

PDA ROLL DISPENSER – 42 patients

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 42 bins H70 L.1/6 + roll dispenser + divider
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



| Reference | Description |
|-------------------|------------------------------------------------------------------------------------------------|
| 8YM64H100DEV42_RM | NEOP "Premium" medicine trolley 600x400 H100 PDA roll dispenser H70 42 patients - magenta pink |
| 8YM64H100DEV42_BC | NEOP "Premium" medicine trolley 600x400 H100 PDA roll dispenser H70 42 patients - sky blue |
| 8YM64H100DEV42_VP | NEOP "Premium" medicine trolley 600x400 H100 PDA roll dispenser H70 42 patients - apple green |


Label covers to order if necessary

Trolleys height 120 cm

PDA ROLL DISPENSER – 36 patients

600 x 400 medicine trolley H120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 36 bins H110 L.1/6 + roll dispenser + divider
- 7 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket




| Reference | Description |
|-------------------|-------------------------------------------------------------------------------------------------|
| 8YM64H120DEV36_RM | NEOP "Premium" medicine trolley 600x400 H120 PDA roll dispenser H110 36 patients - magenta pink |
| 8YM64H120DEV36_BC | NEOP "Premium" medicine trolley 600x400 H120 PDA roll dispenser H110 36 patients - sky blue |
| 8YM64H120DEV36_VP | NEOP "Premium" medicine trolley 600x400 H120 PDA roll dispenser H110 36 patients - apple green |

Label covers to order if necessary

PDA ROLL DISPENSER – 42 patients

600 x 400 medicine trolley H120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 42 bins H110 L.1/6 roll dispenser + divider
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket




| Reference | Description |
|-------------------|-------------------------------------------------------------------------------------------------|
| 8YM64H120DEV42_RM | NEOP "Premium" medicine trolley 600x400 H120 PDA roll dispenser H110 42 patients – magenta pink |
| 8YM64H120DEV42_BC | NEOP "Premium" medicine trolley 600x400 H120 PDA roll dispenser H110 42 patients - sky blue |
| 8YM64H120DEV42_VP | NEOP "Premium" medicine trolley 600x400 H120 PDA roll dispenser H110 42 patients - apple green |

Label covers to order if necessary

PDA ROLL DISPENSER – 54 patients

600 x 400 medicine trolley H120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 54 bins H70 L.1/6 roll dispenser + divider
- 10 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H115
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



| Reference | Description |
|-------------------|------------------------------------------------------------------------------------------------|
| 8YM64H120DEV54_RM | NEOP "Premium" medicine trolley 600x400 H120 PDA roll dispenser H70 54 patients – magenta pink |
| 8YM64H120DEV54_BC | NEOP "Premium" medicine trolley 600x400 H120 PDA roll dispenser H70 54 patients - sky blue |
| 8YM64H120DEV54_VP | NEOP "Premium" medicine trolley 600x400 H120 PDA roll dispenser H70 54 patients - apple green |

Label covers to order if necessary

Standard configurations – BASIC

Structure WITHOUT internal fittings or peripheral accessories.

Trolleys height 100 cm


TROLLEY WITHOUT FOLDABLE EXTENSION SHELF

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 3 standard rails (1 top right, 1 bottom right, 1 bottom left)

Trolley with standard dividers - 40 mm pitch

| Reference | Description |
|--------------------|----------------------------------------------------------------------------------------------|
| 8YCHM164H100FRC_RM | NEOP "Basic" medicine trolley 600x400 H100 Roller shutter door + code lock – magenta pink |
| 8YCHM164H100FRC_BC | NEOP "Basic" medicine trolley 600x400 H100 Roller shutter door + code lock – sky blue |
| 8YCHM164H100FRC_VP | NEOP "Basic" medicine trolley 600x400 H100 Roller shutter door + code lock – apple green |



Trolley with space-optimized dividers* – 62 mm pitch

| Reference | Description |
|-----------------------|----------------------------------------------------------------------------------|
| 8YCHM164H100FRC_RM_ED | NEOP "Basic" medicine trolley 600x400 H100 P62 Roller + code lock – magenta pink |
| 8YCHM164H100FRC_BC_ED | NEOP "Basic" medicine trolley 600x400 H100 P62 Roller + code lock – sky blue |
| 8YCHM164H100FRC_VP_ED | NEOP "Basic" medicine trolley 600x400 H100 P62 Roller + code lock – apple green |


TROLLEY WITH FOLDABLE EXTENSION SHELF

600 x 400 medicine trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 standard rails (1 top right, 1 bottom right, 1 bottom left)
- 1 foldable extension shelf

Trolley with standard dividers - 40 mm pitch

| Reference | Description |
|--------------------|----------------------------------------------------------------------------------------------------------------------------|
| 8YCHM164H100RCT_RM | NEOP "Basic" medicine trolley 600x400 H100 roller shutter door + code lock + foldable extension shelf – magenta pink |
| 8YCHM164H100RCT_BC | NEOP "Basic" medicine trolley 600x400 H100 roller shutter door + code lock + foldable extension shelf – sky blue |
| 8YCHM164H100RCT_VP | NEOP "Basic" medicine trolley 600x400 H100 roller shutter door + code lock + foldable extension shelf – apple green |



Trolley with space-optimized dividers* – 62 mm pitch

| Reference | Description |
|-----------------------|-------------------------------------------------------------------------------------------------------|
| 8YCHM164H100FRC_RM_ED | NEOP "Basic" medicine trolley 600x400 H100 P62 Roller + code lock + extension shelf – magenta pink |
| 8YCHM164H100FRC_BC_ED | NEOP "Basic" medicine trolley 600x400 H100 P62 Roller + code lock + extension shelf – sky blue |
| 8YCHM164H100FRC_VP_ED | NEOP "Basic" medicine trolley 600x400 H100 P62 Roller + code lock + extension shelf – apple green |

* Optimized for Modulo XL, Modulo pill dispensers Classic 600 x 400 and Eurojour


Trolleys height 120 cm

TROLLEY WITHOUT FOLDABLE EXTENSION SHELF

600 x 400 medicine trolley H120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 3 rails (1 top right, 1 bottom right, 1 bottom left)

| Reference | Description |
|--------------------|-------------------------------------------------------------------------------------------|
| 8YCHM164H120FRC_RM | NEOP "Basic" medicine trolley 600x400 H120 roller shutter door + code lock – magenta pink |
| 8YCHM164H120FRC_BC | NEOP "Basic" medicine trolley 600x400 H120 roller shutter door + code lock – sky blue |
| 8YCHM164H120FRC_VP | NEOP "Basic" medicine trolley 600x400 H120 roller shutter door + code lock – apple green |




TROLLEY WITH FOLDABLE EXTENSION SHELF

600 x 400 medicine trolley H120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 rails (1 bottom right, 1 bottom left)
- 1 foldable extension shelf

| Reference | Description |
|--------------------|----------------------------------------------------------------------------------------------------------------------|
| 8YCHM164H120RCT_RM | NEOP "Basic" medicine trolley 600x400 H120 Roller shutter door + code lock + foldable extension shelf – magenta pink |
| 8YCHM164H120RCT_BC | NEOP "Basic" medicine trolley 600x400 H120 Roller shutter door + code lock + foldable extension shelf – sky blue |
| 8YCHM164H120RCT_VP | NEOP "Basic" medicine trolley 600x400 H120 Roller shutter door + code lock + foldable extension shelf – apple green |



Custom configurations

STEP 1 - Trolley structure

Choice of trolley structure according to inside panel models

| | Trolley height 80 cm* | Trolley height 100 cm* | Trolley height 120 cm* |
|-------------------------------------------------------------------------|-----------------------|------------------------|------------------------|
| Standard inside panels for standard drawers 40 mm level | 8Y164H80P4 | 8Y164H100P4 | 8Y164H120P4 |
| Inside panels for telescopic drawers | - | 8Y164H100T | 8Y164H120T |
| Inside panels for « Adhesia Equipments » pill dispensers 40 mm level | - | 8Y164H100W | 8Y164H120W |

This reference includes:

- Structure
- Wheels
- Inside panels

STEP 2 - Trolley customization

Selection of elements 1 to 4 is compulsory

| | | | Trolley height 80 cm* | Trolley height 100 cm* | Trolley height 120 cm* |
|----------|-----------------------|--------------------------------------|-----------------------|------------------------|------------------------|
| 1 | Upper shelf | thermoplastic | 8Y001_6400TH | 8Y001_6400TH | 8Y001_6400TH |
| 2 | Closing system | Roller shutter door | 8Y002_6008F3 | 8Y002_6010F3 | 8Y002_6012F3 |
| | | Door wing for telescopic drawers | - | 8Y002_0010F2 | 8Y002_0010F2 |
| 3 | Code lock | For roller shutter door | 8Y003CODE3 | 8Y003CODE3 | 8Y003CODE3 |
| | | For door wing for telescopic drawers | - | 8Y003CODE2 | - |
| | Badge lock | For roller shutter door | - | 8Y003BADGE3 | 8Y003BADGE3 |
| | | For door wing for telescopic drawers | - | 8Y003BADGE2 | - |
| 4 | Colour | Magenta pink | 8Y004_0008RM | 8Y004_0010RM | 8Y004_0012RM |
| | | Apple green | 8Y004_0008VP | 8Y004_0010VP | 8Y004_0012VP |
| | | Sky blue | 8Y004_0008BC | 8Y004_0010BC | 8Y004_0012BC |

*For precise height of working surface, see technical data



Thermoplastic upper shelf with ergonomic gripping area



Roller shutter closure



Door wing



Code lock



Badge lock

Magenta pink

Apple green

Sky blue


Accessories and internal fittings

Extension shelf, standard rail, ergonomic handle, motorization, arm for laptop, upper frame, drawer...

Please refer to product data sheets.



Technical data



Materials

| | |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Structure | Complete welded one-piece structure, EPOXY anti-corrosion treated H 80: Steel (electro-galvanized) / H 100: Aluminium / H 120: Aluminium |
| Work surface | Thermoplastic: ABS / PMMA |
| Slides | ABS (Acrylonitrile Butadiene Styrene) or ballbearing slide made of galvanized steel with shock absorber |
| Closing system | Roller shutter door: PP (Polypropylene), 2 mm thickness Door wing: Aluminium with EPOXY painting |
| Wheels | Simple non marking rubber wheels - Polyamide cover (PA) |
| Bumpers | 4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate) |
| Strip of colour | 70 micron calendered PVC vinyl film, with acrylic adhesive |

Additional information

| | | | |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Wheels | 125 mm diameter - 2 breaking wheels + 2 swivelling wheels - wheel and wheel-pin double blocking | | |
| Empty weight | H 80: 27 kg | H 100: 24 kg | H 120: 26 kg |
| Capacity | H 80 with standard dividers (40mm): - 15 levels of 40 mm for trolley with roller shutter door - 16 levels of 40 mm for trolley without roller shutter door | H 100 with standard dividers (40mm): - 19 levels for trolley with roller shutter - 20 levels without roller shutter door H 100 with dividers at optimized size (62 mm): - 11 levels without roller shutter door - 12 levels without roller shutter door H 100 with walls for drawers telescopic: - 10 telescopic drawers H 80 mm | H 120 with standard dividers (40mm): - 25 levels for trolley with roller - 26 levels for trolley without roller shutter door H 120 with walls for drawers telescopic: - 13 telescopic drawers H 80 mm |
| Locking system | Code lock / Badge lock | | |

Packaging

| | | | |
|-------------------------------------------|----------------------------------------------------------------------------------------|-----------------------------|---------------------------------------------------------------------------------|
| Packing | Individually on pallet strapped and film-wrapped, or on wheels, with proper protection | Dimensions L X D x H | H 80: 80 x 120 x 120 cm H 100: 80 x 120 x 145 cm H 120: 80 x 120 x 160 cm |
| Gross weight (trolley + packaging) | H 80: 39 kg H 100: 36 kg H 120: 38 kg | Volume | 1,15 m ³ |