

# Neop | 600 x 400 emergency trolleys with joint drawers

## Material and products immediately available for first aid

- Robustness: excellent resistance, one-piece structure made of aluminium
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ward name identification
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours



Robust



Hygienic



Ergonomic



Customizable

- Height 100 cm
- Ergonomic gripping area
- Door wing lock with 1 batch of 20 seals. Access to the content of the trolley breaking the seal by simple twist. After each use, the trolley is reequipped, then closed by a single-use numbered seal.
- 5 joint drawers
- Cardiac massage board on support

### Colour

Garance red



# Standard configuration - PREMIUM

600 x 400 emergency trolleys with joint drawers

Ref. 8Y164H100P4.U2


## Standard equipments

### Peripheral accessories

- 1 batch of 20 seals
- 1 cardiac massage board holder
- 1 cardiac massage board
- 1 monitoring holder
- 1 IV pole with 2 hooks
- 1 lateral standard rail - bottom right
- 1 oxygen tank holder

### Internal fittings

- 1 drawer 600 x 400 height 75 mm
- 3 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm



## Accessories and internal fittings

Ergonomic handle, upper frame, standard rail, bag holder...  
*Please refer to product data sheets.*



# Technical data

## Overall dimensions of the trolley



Les cotes hors tout sont relevées avec les roues dans la position la moins encombrante.

## Useful dimensions of the worktop



## Worktop height

H. 1057 mm

## Materials

Structure	Complete welded and EPOXY anti-corrosion treated one-piece structure
Work surface	ABS / PMMA thermoplastic
Closing system	Door wing: Aluminium with EPOXY painting
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)
Strip of colour	70 micron calendered PVC vinyl film with acrylic adhesive

## Additional information

Wheels	125 mm diameter - 2 breaking wheels + 2 swivelling wheels - wheel and wheel-pin double blocking
Empty weight	24 kg
Capacity	20 levels of 40 mm
Locking system	Seals

## Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	Dimensions L X D x H	L. 80 x D. 120 x H. 145 cm
Gross Weight (trolley + packaging)	36 kg	Volume	1,15 m <sup>3</sup>