

Neop | 600 x 400 nursing trolleys

Care management

- Robustness: excellent resistance, one-piece structure made of steel (height 80 cm) and of aluminium (height 100 cm)
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ergonomic and handy: designed for medical staff comfort
- Ward name identification



Robust



Hygienic



Ergonomic



Customizable

- Single-column trolley
- 2 heights: 80 cm or 100 cm
- Roller shutter door or door wing
- Code lock or badge lock
- 2 available inside panel models (slides) according to internal fittings

Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue



Standard configurations - PREMIUM


Structure WITH internal fittings and peripheral accessories.
Remember to add the color reference of the trolley.

JOINT DRAWERS WITHOUT CLOSURE

1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- 1 drawer 600 x 400 height 75 mm
- 3 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm
- 2 rails (down right and down left)



Reference

8YS64H100TJ7

JOINT DRAWERS WITHOUT CLOSURE + ACC.

1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- 1 drawer 600 x 400 height 75 mm
- 3 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm
- 2 rails (down right et middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin



Reference

8YS64H100TJ4SF

JOINT DRAWERS WITHOUT CLOSURE + 2 RAILS

1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 1 drawer 600 x 400 height 75 mm
- 3 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm
- 2 rails (down right et middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin



Reference

8YS64H100TJ4

JOINT DRAWERS WITH CLOSURE + 3 RAILS

1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 1 drawer 600 x 400 height 75 mm
- 3 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm
- 3 rails (top right, down right, down left)
- 1 standard support for sharp container
- 1 double bag holder
- 1 large side basket



JOINT DRAWERS WITHOUT CLOSURE

3 small / 2 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- 3 drawers 600 x 400 height 75 mm
- 2 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm
- 2 rails (down right et middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin




JOINT DRAWERS WITH CLOSURE + 2 RAILS

3 small / 2 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 drawers 600 x 400 height 75 mm
- 2 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm
- 2 rails (down right et middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin



JOINT DRAWERS WITH CLOSURE + 3 RAILS 3 small / 2 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 drawers 600 x 400 height 75 mm
- 2 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm
- 3 rails (top right, down left, middle left)
- 1 standard support for sharp container
- 1 double bag holder
- 1 large side basket



Reference

8YS64H100TJ2

JOINT DRAWERS WITH CLOSURE + 2 RAILS 4 medium / 2 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 4 drawers 600 x 400 height 115 mm
- 2 drawers 600 x 400 height 155 mm
- 2 rails (top right and et middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin



Reference

8YS64H100TJ6

JOINT DRAWERS WITH CLOSURE + 3 RAILS 4 medium / 2 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 4 drawers 600 x 400 height 115 mm
- 2 drawers 600 x 400 height 155 mm
- 3 rails (top right, down right, down left)
- 1 standard support for sharp container
- 1 double bag holder
- 1 large side basket



Reference

8YS64H100TJ3

TELESCOPIC DRAWERS WITHOUT CLOSURE 5 medium

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- 5 telescopic drawers H. 160 mm with ISO tray H. 100 mm
- 2 rails (down right and middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin




Reference

8YS64H100TT6SF

TELESCOPIC DRAWERS WITH CLOSURE + 2 RAILS 5 medium

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 5 telescopic drawers H. 160 mm with ISO tray H. 100 mm
- 2 rails (down right and middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin



Reference

8YS64H100TT6

TELESCOPIC DRAWERS WITH CLOSURE + 3 RAILS 5 medium

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 5 telescopic drawers H. 160 mm with ISO tray H. 100 mm
- 3 rails (top right, down right, down left)
- 1 standard support for sharp container
- 1 double bag holder
- 1 large side basket



Reference

8YS64H100TT1

TELESCOPIC DRAWERS WITHOUT CLOSURE

1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- 1 telescopic drawer H. 80 mm with ISO tray H. 50 mm
- 3 telescopic drawers H. 160 mm with ISO tray H. 100 mm
- 1 telescopic drawer H. 240 mm with ISO tray H. 200 mm



Reference


8YS64H100TT7

TELESCOPIC DRAWERS WITHOUT CLOSURE + 2 RAILS

1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- 1 telescopic drawer H. 80 mm with ISO tray H. 50 mm
- 3 telescopic drawers H. 160 mm with ISO tray H. 100 mm
- 1 telescopic drawer H. 240 mm with ISO tray H. 200 mm
- 2 rails (down right and middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin



Reference

8YS64H100TT4SF

TELESCOPIC DRAWERS WITH CLOSURE + 2 RAILS

1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 1 telescopic drawer H. 80 mm with ISO tray H. 50 mm
- 3 telescopic drawers H. 160 mm with ISO tray H. 100 mm
- 1 telescopic drawer H. 240 mm with ISO tray H. 200 mm
- 2 rails (down right and middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin



Reference

8YS64H100TT4

TELESCOPIC DRAWERS WITH CLOSURE + 3 RAILS

1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 1 telescopic drawer H. 80 mm with ISO tray H. 50 mm
- 3 telescopic drawers H. 160 mm with ISO tray H. 100 mm
- 1 telescopic drawer H. 240 mm with ISO tray H. 200 mm
- 3 rails (top right, down right, down left)
- 1 standard support for sharp container
- 1 double bag holder
- 1 large side basket



Reference

8YS64H100TT2

TELESCOPIC DRAWERS WITH CLOSURE + 2 RAILS

3 small / 2 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 telescopic drawers H. 80 mm with ISO tray H. 50 mm
- 2 telescopic drawers H. 160 mm with ISO tray H. 100 mm
- 1 telescopic drawer H. 240 mm with ISO tray H. 200 mm
- 2 rails (down right and middle left)
- 1 standard support for sharp container
- 1 simple mobile dustbin



Reference


8YS64H100TT5

TELESCOPIC DRAWERS WITH CLOSURE + 3 RAILS

3 small / 2 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 telescopic drawers H. 80 mm with ISO tray H. 50 mm
- 2 telescopic drawers H. 160 mm with ISO tray H. 100 mm
- 1 telescopic drawer H. 240 mm with ISO tray H. 200 mm
- 3 rails (top right, down right, down left)
- 1 standard support for sharp container
- 1 double bag holder
- 1 large side basket



Reference

8YS64H100TT3

Standard configurations – BASIC

Structure WITHOUT internal fittings or peripheral accessories.
Remember to add the color reference of the trolley.

Trolleys height 80 cm


Reference	Description
8YCHM164H80FBC	NEOP Basic 600x400 H80 nursing trolley with thermoplastic upper shelf + door wing + code lock + standard inside panels (40 mm level)
8YCHM164H80CFBC	NEOP Basic 600x400 H80 nursing trolley with compact upper shelf + door wing + code lock + standard inside panels (40 mm level)

Trolleys height 100 cm

TROLLEY WITH COMPACT UPPER SHELF

600 x 400 nursing trolley H100 equipped with:

- Compact upper shelf
- Door wing
- Code lock




Reference	Description
8YCHM164H100CFBC	NEOP Basic 600x400 H100 nursing trolley with compact upper shelf + door wing + code lock + standard inside panels (40 mm level)

TROLLEY WITHOUT FOLDABLE EXTENSION SHELF

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 rails (top right, down right and down left)



Reference	Description
8YCHS164H100	NEOP Basic 600x400 H100 nursing trolley + door wing + code lock standard inside panels (40 mm level)
8YCHST164H100	NEOP Basic 600x400 H100 nursing trolley + door wing + code lock inside panels for telescopic drawers

TROLLEY WITH FOLDABLE EXTENSION SHELF

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 2 rails (down right and down left)
- 1 foldable extension shelf




Reference	Description
8YCHS164H100T	NEOP Basic 600x400 H100 nursing trolley + door wing + code lock + foldable extension shelf + standard inside panels (40 mm level)
8YCHST164H100T	NEOP Basic 600x400 H100 nursing trolley + door wing + code lock + foldable extension shelf + inside panels for telescopic drawers

TROLLEY WITH ERGONOMIC HANDLE

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- Standard inside panels (40 mm level)
- 1 ergonomic handle
- 3 rails (middle left, down left and down right)




Reference	Description
8YCHM164H100_BCP	NEOP Basic 600x400 H100 nursing trolley + door wing + code lock + 3 rails + ergonomic handle + standard inside panels (40 mm level)

TROLLEY WITH MOTORIZATION OPTION

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- Standard inside panels (40 mm level)
- 3 rails (middle left, down left and down right)
- Motorization system



Reference	Description
8YCHM164H100_BCM	NEOP Basic 600x400 H100 nursing trolley with motorization option + door wing + code lock + 3 rails + standard inside panels (40 mm level)

Custom configurations

STEP 1 - Trolley structure

Choice of trolley structure according to inside panel model.

	Trolley height 80 cm*	Trolley height 100 cm*	Trolley height 120 cm*
Standard inside panels (40 mm level)	8Y164H80P4	8Y164H100P4	8Y164H120P4
Inside panels for telescopic drawers	-	8Y164H100T	8Y164H120T

This reference includes the structure, the wheels and the inside panels.

STEP 2 - Trolley customization

		Trolley height 80 cm*	Trolley height 100 cm*	Trolley height 120 cm*	
1	Upper shelf	Thermoplastic with lock	8Y001_6400TH		
		Thermoplastic without lock	8Y001_6400THSP		
		Compact laminate	8Y001_6400CP1		
2	Closing system	Roller shutter door for standard inside panels	8Y002_6008F3	8Y002_6010F3	8Y002_6012F3
		Door wing for telescopic drawers	-	8Y002_0010F2	8Y002_0010F2
		Door wing for standard inside panels	8Y002_0008F1	8Y002_0010F1	-
3	Code lock	For roller shutter door	8Y003CODE3		
		For door wing for telescopic drawers	8Y003CODE2		
		For door wing for standard inside panels	8Y003CODE1		
	Badge lock	For roller shutter door	8Y003BADGE3		
		For door wing for telescopic drawers	8Y003BADGE2		
		For door wing for standard inside panels	8Y003BADGE1		
4	Colour	Magenta pink	8Y004_0008RM	8Y004_0010RM	8Y004_0012RM
		Apple green	8Y004_0008VP	8Y004_0010VP	8Y004_0012VP
		Sky blue	8Y004_0008BC	8Y004_0010BC	8Y004_0012BC
5	Marking	Specific marking			8Y132



Thermoplastic upper shelf with ergonomic gripping area



Roller shutter closure



Code lock

Magenta pink



Upper shelf made of compact laminate



Door wing



Badge lock

Apple green

Sky blue

Accessories and internal fittings

Extension shelf, ergonomic handle, , arm for laptop, upper frame, dustbin, drawer, motorization...



*For precise height of working surface, see technical data

Technical data

Trolley with thermoformed plastic top shelf

Overall dimensions of the trolley



Useful dimensions of the worktop



Worktop height

H 120 : 1297 mm
H 100 : 1057 mm
H 80 : 897 mm

Trolley with top shelf in compact laminate

Overall dimensions of the trolley



Useful dimensions of the worktop



Worktop height

H 120 : 1304 mm
H 100 : 1064 mm
H 80 : 904 mm

The overall dimensions are measured with the wheels in the least bulky position.
The depth incorporates a lock by default.

Materials

Structure	Complete welded one piece structure, EPOXY anti corrosion treated (H 80: steel electro galvanized) / H 100: aluminium / H 120: aluminium
Work Surface	Thermo plastic: ABS / PMMA or compact laminate 10 mm thickness
Slides	ABS (Acrylonitrile Butadiene Styrene) or ballbearing slides made of galvanized steel with shock absorber
Closing System	Roller shutter door : PP Polypropylene), 2 mm thickness Door wing : Aluminium with EPOXY painting
Wheels	Simple non marking rubber wheels Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)
Strip of colour	70 micron calendered PVC vinyl film, with acrylic adhesive

Additional information

Wheels	125 mm diameter 2 breaking wheels + 2 swivelling wheels wheel and wheel pin double blocking		
Empty weight	H 80: 27 kg	H 100: 24 kg	H 120: 26 kg
Capacity	H 80 with standard inside panels (40 mm level): - 15 levels for trolley with roller shutter door - 16 levels for trolley without roller shutter door or trolley with door wing	H 100 with standard inside panels (40 mm level): - 19 levels for trolley with roller shutter door - 20 levels for trolley without roller shutter door or trolley with door wing H 100 with optimized inside panels (62 mm level): - 11 levels for trolley with roller shutter door - 12 levels for trolley without roller shutter door H 100 with inside panels for telescopic drawers: 10 drawers H 80 mm	H 120 with standard inside panels (40 mm level): - 25 levels for trolley with roller shutter door - 26 levels for trolley without roller shutter door H 120 with inside panels for telescopic drawers: 13 drawers H 80 mm
Locking system	Code lock (with spare key) or badge lock		

Packaging

Packing	Individually on pallet strapped and film wrapped, or on wheels, with proper protection
Gross weight (trolley + packaging)	H 80: 39 kg H 100: 36 kg H 120: 38 kg
Dimensions	H 80: L 0,80 x D 1,20 x H 1,20 m H 100: L 0,80 x D 1,20 x H 1,45 m H 120: L 0,80 x D 1,20 x H 1,60 m
Volume	1,15 m ³