

## An optimal configuration

- Robustness: excellent resistance, one-piece structure made of aluminium
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ward name identification
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours



Robust



Hygienic



Ergonomic



Customizable

- Single-column trolley 600x400
- Height: 100 cm
- Ergonomic gripping area

### Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue




# Standard configuration - PREMIUM

## Anesthesia trolley 600 x 400 H. 100 mm

Trolley 600 x 400 H100 equipped with:

- Thermoplastic upper shelf
- 3 drawers 600 x 400 H. 75 mm with compartment kit 3x2
- 2 drawers 600 x 400 H. 155 mm with compartment kit 2x2
- 1 drawer 600 x 400 H. 235 mm
- 1 horizontal handle
- 2 lateral rails (top right - bottom right)
- 1 adjustable sharp disposal container holder
- 1 catheter holder
- 1 upper frame equipped with:
  - 2 lateral rails
  - 6 tipping container L. 601 x D. 93 x H. 114 mm
  - 5 tipping container L. 601 x D. 133 x H. 166 mm



Reference	Description
8YCHANSS164H100_RM	Neop anesthesia trolley 600x400 H100 - magenta pink
8YCHANSS164H100_VP	Neop anesthesia trolley 600x400 H100 - apple green
8YCHANSS164H100_BC	Neop anesthesia trolley 600x400 H100 - sky blue

*\*For precise height of working surface, see next page*

## Accessories

Extension shelf, standard rail, arm for laptop, bag holder...

*Please refer to product data sheets.*



# Technical data

Overall dimensions of the trolley



Useful dimensions of the worktop



Worktop height

1057 mm

## Materials

Structure	Complete welded one-piece aluminium structure, EPOXY anti-corrosion treated
Work surface	Thermoplastic: ABS / PMMA
Slides	ABS (Acrylonitrile Butadiene Styrene)
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)
Strip of colours	70 micron calendered PVC vinyl film, with acrylic adhesive

## Additional information

Wheels	Diameter 125 mm - 2 breaking wheels + 2 swivelling wheels Wheel and wheel-pin double blocking
Empty weight	24 kg
Capacity	20 levels of 40 mm

## Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	Dimensions	L 0,80 x D 1,20 x H 1,45 m
Gross weight (trolley + packaging)	36 kg	Volume	1,15 m <sup>3</sup>