


SATELLITE

MULTI-SENSORY NURSING TROLLEY



PRATICIMA
healthcare solutions



▶ Video

SATELLITE

Unique design for multi-sensory care

- A cart which combines all kind of multi-sensory elements
- A mobile solution which doesn't need any dedicated room
- Usable in any facility including wet environment: bathroom or balneotherapy
- Adapted to people who have difficulty moving: bedridden or in (wheel) chair

An entire multi-sensory experience based on:



HEARING



SIGHT



SMELL



TOUCH

Musical speakers
high definition


Adjustable spotlight
for video broadcasting
on wall or ceiling



Essential oil diffuser:
6 different fragrances

Fiber optic cables,
with colour changes,
located at optimal
height for people in
bed or in (wheel) chair

**Storage drawers
of small materials:**
light balls, hedgehog balls,
heating pad...

Water bubble panel
big size: seamless flow
and quiet pump





Non-drug therapy as a complement to healthcare pathway


Its benefits

- It prevents and fights against psycho-behavioral disorders
- It provides time and space benchmarks
- It promotes relaxation and well-being
- It helps fight against pain
- It appeals to sensory and emotional memory
- It enables to enjoy customized care adapted to the moment and state of the users, thanks to **6 predefined atmospheres or a free program**

SATELLITE

A compact trolley, ideal for facilities
in search of space saving

DESIGNED IN COOPERATION WITH DOCTORS AND OCCUPATIONAL THERAPISTS



SATELLITE

A revolution to provide multi-sensory care

An easy-to-use ergonomic trolley

- Care programmes can start by a single action: from **control panel** or remote control, from **touchpad**
- **Care implementation is simple and quick**, for a higher availability of medical staff
- Aluminium **solid** construction
- **Autonomous**, equipped with batteries
- **Compact**, offers lightweight handling

6 PRE-SET CARE PROGRAMMES



STIMULATING



BALNEO



BEFORE MEAL



SOOTHING



RELAXING



NIGHT

Simultaneous activation of a video, a music, a fragrance, the water bubble panel and of the fiber optic cables

CARE CUSTOMIZATION WITH THE « FREE ATMOSPHERE » FUNCTION

Selection of the music, video, fragrance, with or without activation of the fiber optic cables and water bubble panel



PRATICIMA
healthcare solutions

Rue du Pou du Ciel
ZI de Reyrieux – CS20129 Reyrieux
01600 TRÉVOUX – FRANCE

Tel: +33 (0)4 74 00 98 11
Fax: +33 (0)4 74 00 98 36
E-mail: praticima@praticima.fr

www.praticima.fr

