


Medicine dispensing Pil 7[®] dispensing trolley

The Pil 7[®] range of medicine trolleys are ideal for distribution in institutions (elderly homes, nursing homes...).

Our **pill dispensers with removable dividers** are totally adapted for medication for autonomous people. Their ergonomic design is an advantage for daily use.



Hygienic



Ergonomic

Key points

- ▶ Storage, transfer and distribution in one trolley only
- ▶ Full opening of the dispenser for medicine filling
- ▶ Resistant to the disinfecting products of hospital environments

Medicine dispensing

Pil 7[®] dispensing trolley

Pil 7[®] dispensing trolley

- ▶ Dimensions: L. 615 x D. 468 x H. 1128 mm
- ▶ Capacity: 30 piluliers
- ▶ Anodised aluminium structure
- ▶ ABS side walls, coloured body
- ▶ White thermoplastic work surface with rounded raised edges
- ▶ Curtain door with tool-free removable slats
- ▶ Locking system, code lock
- ▶ 4 polyamide swivelling wheels Ø 100 mm, 1 with brake
- ▶ 1 handle
- ▶ 4 standard rails

Description

Reference

Pil 7[®] dispensing trolley

8D030_PIL7_CODE



Pil 7[®] Evolution pill dispenser

- ▶ Weekly pill dispenser of 28 compartments which enables to take medicines 4 times a day during one week
- ▶ Mobile and adjustable dividers according to prescription
- ▶ 3 different colours for dividers which enable a better identification
- ▶ The 7 days of the week are printed
- ▶ Mentions: « Morning - Midday - Evening - Night » on dividers
- ▶ Transparent sliding lids for a rapid use
- ▶ Full opening of the dispenser for medicine filling
- ▶ Easy access into each compartment
- ▶ Dimensions: L. 320 x D. 200 x H. 40 mm

Description

Packing

Reference

Evolution Pil 7[®] pill dispenser without end stops

Unit

CD23605

Evolution Pil 7[®] pill dispenser without end stops

Batch of 10

CD23606

