

Neop | 400 x 400 single-column nursing trolleys

Care management system

- Robustness: excellent resistance, one-piece structure made of aluminium
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ward name identification
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours



Robust



Hygienic



Ergonomic



Customizable

- Single-column trolley
- Height* 100 cm
- Roller shutter door or door wing
- Code lock (with spare key) or badge lock
- 2 available inside panel models (slides) according to required internal fittings
- Ergonomic gripping area

Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue



*For precise height of working surface, see « Technical data »

Standard configurations - PREMIUM

Structure WITH internal fittings and peripheral accessories.


Rails position:

- top right
- bottom right
- bottom left

TELESCOPIC DRAWERS – 5 medium

400 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 5 telescopic drawers F400 H160 with perforated ISO tray H100
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket




Reference	Description
8YS44H100TT1_RM	NEOP "Premium" 400x400 nursing trolley H100 telescopic drawers 5M – magenta pink
8YS44H100TT1_BC	NEOP "Premium" 400x400 nursing trolley H100 telescopic drawers 5M – sky blue
8YS44H100TT1_VP	NEOP "Premium" 400x400 nursing trolley H100 telescopic drawers 5M – apple green

TELESCOPIC DRAWERS – 1 small / 3 medium / 1 large

400 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 1 telescopic drawer F400 H80 with standard ISO tray H50
- 3 telescopic drawers F400 H160 with perforated ISO tray H100
- 1 telescopic drawer F400 H240 with perforated ISO tray H100
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket




Reference	Description
8YS44H100TT2_RM	NEOP "Premium" 400x400 nursing trolley H100 telescopic drawers 1S / 3M / 1L - magenta pink
8YS44H100TT2_BC	NEOP "Premium" 400x400 nursing trolley H100 telescopic drawers 1S / 3M / 1L - sky blue
8YS44H100TT2_VP	NEOP "Premium" 400x400 nursing trolley H100 telescopic drawers 1S / 3M / 1L - apple green

TELESCOPIC DRAWERS – 3 small / 2 medium / 1 large

400 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 telescopic drawers F400 H80 with standard ISO tray H50
- 2 telescopic drawers F400 H160 with perforated ISO tray H100
- 1 telescopic drawer F400 H240 with perforated ISO tray H100
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



Reference	Description
8YS44H100TT3_RM	NEOP "Premium" 400x400 nursing trolley H100 telescopic drawers 3S / 2M / 1L – magenta pink
8YS44H100TT3_BC	NEOP "Premium" 400x400 nursing trolley H100 telescopic drawers 3S / 2M / 1L – sky blue
8YS44H100TT3_VP	NEOP "Premium" 400x400 nursing trolley H100 telescopic drawers 3S / 2M / 1L – apple green


Standard configurations – BASIC

Structure WITH internal fittings and peripheral accessories.

TROLLEY WITHOUT FOLDABLE EXTENSION SHELF

400 x 400 single-column dressing trolley H100 equipped with:

- Inside panels for telescopic drawers
- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 standard rails (top right, bottom right, bottom left)



Reference	Description
8YCHST144H100_RM	NEOP "Basic" nursing trolley 400x400 H100 door wing + code lock - magenta pink
8YCHST144H100_BC	NEOP "Basic" nursing trolley 400x400 H100 door wing + code lock - sky blue
8YCHST144H100_VP	NEOP "Basic" nursing trolley 400x400 H100 door wing + code lock - apple green

TROLLEY WITH FOLDABLE EXTENSION SHELF

400 x 400 single-column dressing trolley H100 equipped with:

- Inside panels for telescopic drawers
- Thermoplastic upper shelf
- Door wing
- Code lock
- 2 standard rails (bottom right, bottom left)
- 1 foldable extension shelf



Reference	Description
8YCHST144H100T_RM	NEOP "Basic" 400x400 nursing trolley H100 door wing + code lock + foldable extension shelf – magenta pink
8YCHST144H100T_BC	NEOP "Basic" 400x400 nursing trolley H100 door wing + code lock + foldable extension shelf – sky blue
8YCHST144H100T_VP	NEOP "Basic" 400x400 nursing trolley H100 door wing + code lock + foldable extension shelf – apple green

Custom configuration

STEP 1 - Trolley structure

Choice of trolley structure according to inside panel models

Trolley height 100 cm*	
Standard inside panels 40 mm levels	Inside panels for telescopic drawers
8Y144H100P4	8Y144H100T

This reference includes:

- **Structure**
- **Wheels**
- **Inside panels**

STEP 2 - Trolley customization

Selection of elements 1 to 4 is compulsory

Trolley with standard
inside panels
ref. 8Y144H100P4

Trolley with inside panels
for telescopic drawers
ref. 8Y144H100T

1	Upper shelf	thermoplastic	8Y001_4400TH	8Y001_4400TH
2	Closing system	Roller shutter door	8Y002_4010F3	-
		Door wing for telescopic drawers	-	8Y002_0010F2
3	Code lock	For roller shutter door	8Y003CODE3	-
		For door wing for telescopic drawers	-	8Y003CODE2
	Badge lock	For roller shutter door	8Y003BADGE3	-
		For door wing for telescopic drawers	-	8Y003BADGE2
4	Colour	Magenta pink	8Y004_0010RM	8Y004_0010RM
		Apple green	8Y004_0010VP	8Y004_0010VP
		Sky blue	8Y004_0010BC	8Y004_0010BC

*For precise height of working surface, see next page



Thermoplastic upper shelf
with ergonomic gripping area



Roller shutter closure



Code lock

Magenta pink

Apple green



Door wing



Badge lock

Sky blue

Accessories and internal fittings

Extension shelf, standard rail, ergonomic handle, motorization, arm for laptop, upper frame, drawer...
Please refer to product data sheets.



Technical data

Overall dimensions of the trolley



Useful dimensions of the worktop



Thermoplastic upper shelf
with ergonomic gripping area

Worktop height

H. 1057 mm

Materials

Structure	Complete welded one-piece aluminium structure, EPOXY anti-corrosion treated
Work surface	Thermoplastic: ABS / PMMA
Slides	ABS (Acrylonitrile Butadiene Styrene) or ballbearing slide made of galvanized steel with shock absorber
Closing system	Roller shutter door: PP (Polypropylene), 2 mm thickness Door wing: Aluminium with EPOXY painting
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)
Strip of colour	70 micron calendered PVC vinyl film, with acrylic adhesive

Additional information

Wheels	125 mm diameter - 2 breaking wheels + 2 swivelling wheels Wheel and wheel-pin double blocking
Empty weight	20 kg
Capacity	- 19 levels of 40 mm for trolley with roller shutter door - 20 levels of 40 mm for trolley without roller shutter door - 10 telescopic drawers H 80 mm
Locking system	Code lock / Badge lock

Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	Dimensions L X D x H	80 x 120 x 120 cm
Gross weight (trolley + packaging)	32 kg	Volume	1,15 m ³