

Neop | 600 x 400 plaster cart

Storage of all instruments and materials which are necessary to make a plaster

- Robustness: excellent resistance, one-piece structure made of aluminium
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ward name identification
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours



Robust



Hygienic



Ergonomic



Customizable

- Height: 100 cm
- Ergonomic gripping area
- Equipped with five drawers with full extension slides and with a lateral support soaking tank
- Many other accessories are available, particularly, the upper frame with tubular jersey roll support



Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue

Cart configuration

STEP 1 - Cart structure

Trolley height 100 cm*
8Y164H100T.PL

This reference includes:

- 5 drawers with full extension slides (1 small, 3 medium and 1 big)
- 1 lateral support for 19 liter soaking tank

STEP 2 - Colour choice

Colour
Magenta pink: 8Y004_0010RM Apple green: 8Y004_0010VP Sky blue: 8Y004_0010BC

Magenta pink

Apple green

Sky blue

STEP 3 - Selection of accessories

Description	Reference
Foldable shelf	8Y008_0400
Upper frame with jersey roll support Made of powder coated steel sheet, colour white Sliding door made of transparent plastic Equipped with 3 dividers mounted on adhesive guides Dimensions: L. 606 x D. 310 x H. 320 mm	8Y012

Standard rail, simple bag holder, double bag holder, ergonomic handle...

Please refer to product data sheets.



Technical data

Overall dimensions of the trolley

H. max : 1830 mm
H. min : 1380 mm

Les cotes hors tout sont relevées avec les roues dans la position la moins encombrante.



Useful dimensions of the worktop



Tablette supérieure en plastique thermoformé avec zone de préhension ergonomique

Worktop height

H. 1057 mm

Materials

Structure	Complete welded and EPOXY anti-corrosion treated aluminium one-piece structure
Work surface	ABS / PMMA thermoplastic
Slides	Ballbearing slides made of galvanized steel with shock absorber
Closing system	Drawers
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)

Additional Information

Wheels	125 mm diameter - 2 breaking wheels + 2 swivelling wheels Wheel and wheel-pin double blocking
Empty weight	27 kg
Capacity	5 compact drawers: 1 big, 3 medium and 1 small
Locking system	None
Strip of colour	70 micron calendered PVC vinyl film with acrylic adhesive
Equipment	Bowl holder + bowl

Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	Dimensions L X D x H	80 x 120 x 174 cm
Gross weight (trolley + packaging)	39 kg	Volume	1,67 m ³