

Care furniture

600 x 400 Resine[®] trolleys

Our composite trolleys have been designed to satisfy special requirements of health establishments. They are completely modular thanks to a **wide range of accessories**. Their shape and colours make them very beautiful. Our trolleys can be **easily disinfected** and **decontaminated** so they perfectly respect highest safety and hygiene requirements.

Non contractual photos



Hygienic



Ergonomic



Modular

Advantages

- ▶ Ergonomic, handy and light
- ▶ Can be equipped according to needs
- ▶ Resistant to disinfecting products of hospital fields
- ▶ Stain-resistant
- ▶ Entirely recyclable

600 x 400 Resine® trolleys

Technical features

- ▶ Anodised aluminium structure
- ▶ PVC shelves and faces
- ▶ 4 Ø 100 mm wheels, 1 with brake (non-marking rubber)
- ▶ Drawer height: 100 mm
- ▶ Dimensions: L. 700 x D. 500 x H. 985 mm
- ▶ Available colours for drawer handles:

Pop blue	Apple green	Fuschia
12	23	45



2-Shelf trolley


Description	Reference
2-shelf composite trolley + 1 drawer with key lock	FA64201.XX
2-shelf composite trolley + 2 drawers (1 with key lock)	FA64202.XX
2-shelf composite trolley + 3 drawers (1 with key lock)	FA64203.XX

3-Shelf trolley

Description	Reference
3-shelf composite trolley	FA64300.71
3-shelf composite trolley + 1 drawer with key lock	FA64301.XX

Optional equipments

Description	Reference
Drawer divider kit for 600 x 400 drawer	FN03803
3-side bar for 600 x 400 trolley	FA01047



Care furniture

400 x 400 Resine[®] trolleys

Our composite trolleys have been designed to satisfy special requirements of health establishments. They are completely modular thanks to a **wide range of accessories**. Their shape and colours make them very beautiful. Our trolleys can be **easily disinfected** and **decontaminated** so they perfectly respect highest safety and hygiene requirements.



Hygienic



Ergonomic



Modular

Advantages

- ▶ Ergonomic, handy and light
- ▶ Can be equipped according to needs
- ▶ Resistant to disinfecting products of hospital fields
- ▶ Stain-resistant
- ▶ Entirely recyclable

400 x 400 Resine® trolleys

Technical features

- ▶ Anodised aluminium structure
- ▶ PVC shelves and faces
- ▶ 4 Ø 100 mm wheels, 1 with brake (non-marking rubber)
- ▶ Drawer height: 100 mm
- ▶ Dimensions: L. 500 x D. 500 x H. 985 mm
- ▶ Available colours for drawer handles:

Pop blue

12

Apple green

23

Fuschia

45



2- Shelf trolley

Description	Reference
2-shelf composite trolley + 1 drawer with key lock	FA44201.XX
2-shelf composite trolley + 2 drawers (1 with key lock)	FA44202.XX
2-shelf composite trolley + 3 drawers (1 with key lock)	FA44203.XX

3- Shelf trolley

Description	Reference
3-shelf composite trolley	FA44300.71
3-shelf composite trolley + 1 drawer with key lock	FA44301.XX

Optional equipments

Description	Reference
Drawer divider kit for 400 x 400 drawer	FN03829
3-side bar for 400 x 400 trolley	FA01046

