

Neop | 600 x 400 distribution trolleys

Daily distribution of medicines with pill dispensers and nominative trays

- Robustness: excellent resistance, one-piece structure made of steel (height 80 cm) and of aluminium (height 100 and 120 cm)
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Light
- Ergonomic and handy: designed for medical staff comfort
- Modular: wide range of drawers and accessories
- Customizable: ward name identification



Robust



Hygienic



Ergonomic



Customizable

- Single column trolley
- 3 height: 80, 100 or 120 cm
- Roller shutter door or door wing
- Code lock or badge lock
- 4 available inside panel models (slides) according to internal fittings

NEW!

EasyMoove® option: trolley motorization
Neop® Lock option: curtain door motorization

Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue



Standard configurations - PREMIUM

Structure WITH internal fittings and peripheral accessories.
Add the color reference of the trolley.

Trolleys height 100 cm

DIN – 18 patients bins H. 70

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 18 bins H. 70 L.1/3 – white color
- 7 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




Reference

8YM64H100DIN18

DIN – 24 patients bins H. 110

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 24 bins H. 110 L. 1/6 – white color
- 5 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 75 mm
- 1 drawer 600 x 400 H. 115 mm
- 1 rail (bottom right)
- 1 bag holder



Reference

8YM64H100DI24E1

DIN – 24 patients bins H. 70

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 24 bins H. 70 L.1/4 – white color
- 7 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 75
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




Reference

8YM64H100DIN24

DIN – 24 patients bins H. 150

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 1 drawer 600 x 400 H. 75 mm for 60 cups
- 50 transparent cups with lid
- 50 blue cups with lid
- 24 bins L. 1/6 H. 150 mm
- 5 shelves 600 x 400 for bins
- 1 rail (bottom right)
- 1 bag holder



Reference

8YM64H100DI24E2

DIN – 28 patients bins H. 70

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 28 bins H. 70 L. 1/4 – white color
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




Reference

8YM64H100DIN28

PDA PILI – 70 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 batches of 10 PILI dispenser supports
- 7 trays 600 x 400 for PILI delivered without dispenser
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




Reference

8YM64H100PIL70

PDA MODULO CLASSIC – 70 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 7 trays 600 x 400 with 10 Modulo Classic dispenser
- 1 drawer 600 x 400 H115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



Reference

8YM64H100MC70

PDA ROLL DISPENSER – 24 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 24 bins H110 L.1/6 + roll dispenser + divider
- 5 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H75
- 1 drawer 600 x 400 H115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



Reference

8YM64H100DEV24

PDA ROLL DISPENSER – 30 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 30 bins H110 L.1/6 + roll dispenser + divider
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



Reference

8YM64H100DEV30

PDA ROLL DISPENSER – 42 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 42 bins H. 70 L. 1/6 + roll dispenser + divider
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



Reference


8YM64H100DEV42

Trolleys height 120 cm

DIN – 36 patients bins H. 110

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 1 drawer 600 x 400 H. 75 mm for 60 cups
- 50 transparent cups with lid
- 50 blue cups with lid
- 36 bins L. 1/6 H. 110 mm – white color
- 7 tablettes 600 x 400 pour bacs tiroirs
- 1 tiroir 600 x 400 H. 155 mm avec séparateurs
- 1 rail (bottom right)
- 1 bag holder




Reference

8YM64H100DI36E1

PDA ROLL DISPENSER – 36 patients

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 36 bins H. 110 L. 1/6 + roll dispenser + divider
- 7 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




Reference

8YM64H120DEV36

PDA ROLL DISPENSER – 42 patients

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 42 bins H. 110 L. 1/6 roll dispenser + divider
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



Reference

8YM64H120DEV42

PDA ROLL DISPENSER – 54 patients

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 54 bins H. 70 L. 1/6 roll dispenser + divider
- 10 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



Reference

8YM64H120DEV54

Standard configurations – BASIC

Structure WITHOUT internal fittings or peripheral accessories.

Add the color reference of the trolley.

Trolleys height 80 cm


Reference	Description
8YCHM164H80FR	NEOP Basic medicine trolley 600x400 H80 thermoplastic upper shelf + curtain door + standard dividers (40 mm level)
8YCHM164H80CFR	NEOP Basic medicine trolley 600x400 H80 compact upper shelf + curtain door + standard dividers (40 mm level)
8YCHM164H80FRC	NEOP Basic medicine trolley 600x400 H80 thermoplastic upper shelf + curtain door + code lock + standard dividers (40 mm level)
8YCHM164H80CFRC	NEOP Basic medicine trolley 600x400 H80 compact upper shelf + curtain door + code lock + standard dividers (40 mm level)

Trolleys height 100 cm

TROLLEY WITH MANUAL CURTAIN DOOR

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 3 rails (top right, bottom right, bottom left)




Reference	Description
8YCHM164H100FRC	NEOP Basic medicine trolley 600x400 H100 P40 + curtain door + code lock + 3 rails + standard dividers (40 mm level)
8YCHMH100FRC_ED	NEOP Basic medicine trolley 600x400 H100 P62 + curtain door + code lock + 3 rails + space-optimized dividers (62 mm level)

TROLLEY WITH NEOP LOCK® MOTORIZED CURTAIN DOOR

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Closing by motorized curtain
- Badge access management
- Reading lamp
- Bottle holder
- 3 rails (top right, bottom right, top left, bottom left)




Reference	Description
8YCHM164H100RNL	NEOP Basic medicine trolley 600x400 H100 P40 + motorized curtain door + 3 rails + standard dividers (40 mm level)

TROLLEY WITH EASYMOOVE® MOTORIZATION SYSTEM

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 3 rails (middle left, bottom left, bottom right)
- Motorization system




Reference	Description
8YCHM164H100_RCM	NEOP Basic medicine trolley 600x400 H100 P40 + motorization system + curtain door + code lock + 3 rails + standard dividers (40 mm level)
8YCHM164H100P6_RCM	NEOP Basic medicine trolley 600x400 H100 P40 + motorization system + curtain door + code lock + 3 rails + space-optimized dividers (62 mm level)

TROLLEY WITH ERGONOMIC HANDLE

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 1 ergonomic handle
- 3 rails (middle left, bottom right, bottom left)




Reference	Description
8YCHM164H100_RCP	NEOP Basic medicine trolley 600x400 H100 P40 + curtain door + code lock + 3 rails + ergonomic handle + standard dividers (40 mm level)
8YCHM164H100P6_RCP	NEOP Basic medicine trolley 600x400 H100 P62 + curtain door + code lock + 3 rails + ergonomic handle + space-optimized dividers (62 mm level)

TROLLEY WITH SLIDING SHELF

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 rails (bottom right, bottom left)
- 1 sliding extension shelf



Reference	Description
8YCHM164H100RCTC	NEOP Basic medicine trolley 600x400 H100 P40 + curtain door + code lock + sliding shelf + 2 rails + standard dividers (40 mm level)
8YCMH100FRCTC_ED	NEOP Basic medicine trolley 600x400 H100 P62 + curtain door + code lock + sliding shelf + 2 rails + space-optimized dividers (62 mm level)

TROLLEY WITH FOLDABLE SHELF

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 rails (bottom right, bottom left)
- 1 foldable extension shelf



Reference	Description
8YCHM164H100RCT	NEOP Basic medicine trolley 600x400 H100 P40 + curtain door + code lock + foldable shelf + 2 rails + standard dividers (40 mm level)
8YCMH100FRCT_ED	NEOP Basic medicine trolley 600x400 H100 P62 + curtain door + code lock + + foldable shelf + 2 rails + space-optimized dividers (62 mm level)


Trolleys height 120 cm

TROLLEY WITH COMPACT LAMINATE TOP SHELF

600 x 400 medicine trolley H. 120 equipped with:

- Compact laminate top shelf
- Roller shutter door
- Code lock

Reference	Description
8Y164H120P6C_RC	NEOP Basic medicine trolley 600x400 H120 P62 + curtain door + code lock + compact top shelf



TROLLEY WITH THERMOPLASTIC UPPER SHELF

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 3 rails (top right, bottom right, bottom left)

Reference	Description
8YCHM164H120FRC	NEOP Basic medicine trolley 600x400 H120 P40 + curtain door + code lock + 3 rails
8YCHMH120FRC_ED	NEOP Basic medicine trolley 600x400 H120 P62 + curtain door + code lock + 3 rails




TROLLEY WITH FOLDABLE SHELF

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 rails (bottom right, bottom left)
- 1 foldable extension shelf

Reference	Description
8YCHM164H120RCT	NEOP Basic medicine trolley 600x400 H120 P40 + curtain door + code lock + foldable shelf + 2 rails
8YCMH120FRCT_ED	NEOP Basic medicine trolley 600x400 H120 P62 + curtain door + code lock + foldable shelf + 2 rails



Custom configurations

STEP ① - Trolley structure

This reference includes structure, wheels and inside panels.

	Trolley height 80 cm*	Trolley height 100 cm*	Trolley height 120 cm*
Standard inside panels for standard drawers (40 mm level)	8Y164H80P4	8Y164H100P4	8Y164H120P4
Inside panels for telescopic drawers	-	8Y164H100T	8Y164H120T
Inside panels for « Adhesia Equipments » pill dispensers (40 mm level)	-	8Y164H100W	8Y164H120W

STEP ② - Trolley customization

STEP 2 - Trolley customization				Trolley height 80 cm*	Trolley height 100 cm*	Trolley height 120 cm*
1	Upper shelf	Thermoplastic shelf with cutout	8Y001_6400TH			
		Thermoplastic shelf without cutout	8Y001_6400THSP			
		Compact laminate shelf with cutout	8Y001_6400CP1			
2	Closing system	Manual curtain door	8Y002_6008F3	8Y002_6010F3	8Y002_6012F3	
		Motorized curtain door Neop® Lock	-	9U000	-	
		Door wing for telescopic drawers	-	8Y002_0010F2	8Y002_0010F2	
		Door wing for standard dividers	8Y002_0008F1	8Y002_0010F1	-	
3	Code lock	For curtain door	8Y003CODE3			
		For door wing and telescopic drawers	8Y003CODE2			
		For door wing and standard dividers	8Y003CODE1			
	Badge lock	For curtain door	8Y003BADGE3			
		For door wing and telescopic drawers	8Y003BADGE2			
		For door wing and standard dividers	8Y003BADGE1			
4	Colour	Magenta pink	8Y004_0008RM	8Y004_0010RM	8Y004_0012RM	
		Apple green	8Y004_0008VP	8Y004_0010VP	8Y004_0012VP	
		Sky blue	8Y004_0008BC	8Y004_0010BC	8Y004_0012BC	
5	Specific marking		8Y132			
6	EasyMoove® trolley motorization		-	9R002_L01	9R002_L01	

*For precise height of working surface, see technical data

Options, accessories and internal fittings

Neop® Lock curtain door motorization


- 100% motorized curtain door with automatic closing and badge access management
- Avoids repetitive handling: improves working comfort
- Secures the contents of the cart: reduces the risk of theft

EasyMoove® trolley motorization

- 100% motorized driving
- Limits pushing efforts for users
- Reduces the risk of musculoskeletal disorders



Extension shelf, ergonomic handle, IT support, upper frame, drawers...



Technical data

Trolley with thermoformed plastic top shelf

Overall dimensions of the trolley



Useful dimensions of the worktop



Ergonomic gripping area

Worktop height

H. 120 : 1297 mm
H. 100 : 1057 mm
H. 80 : 897 mm

Trolley with compact laminate top shelf

Overall dimensions of the trolley



Useful dimensions of the worktop



Worktop height

H. 120 : 1304 mm
H. 100 : 1064 mm
H. 80 : 904 mm

The overall dimensions are taken with the wheels in the least bulky position.
The depth includes a lock by default.

Materials

Structure	Complete welded one-piece structure, EPOXY anti-corrosion treated H 80: Steel (electro-galvanized) / H 100: Aluminium / H 120: Aluminium
Work surface	Thermoplastic: ABS / PMMA OR compact laminate th.. 10 mm
Slides	ABS (Acrylonitrile Butadiene Styrene) or ballbearing slide made of galvanized steel with shock absorber
Closing system	Roller shutter door: PP (Polypropylene), 2 mm thickness Door wing: Aluminium with EPOXY painting
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)
Strip of colour	70 micron calendered PVC vinyl film, with acrylic adhesive

Additional information

Wheels	125 mm diameter - 2 breaking wheels + 2 swivelling wheels - wheel and wheel-pin double blocking		
Empty weight	H. 80 : 27 kg	H. 100 : 24 kg	H. 120 : 26 kg
Capacity	H 80 with standard dividers (40mm): - 15 levels of 40 mm for trolley with roller shutter door - 16 levels of 40 mm for trolley without roller shutter door	H 100 with standard dividers (40 mm): - 19 levels for trolley with roller shutter - 20 levels without roller shutter door H 100 with dividers at optimized size (62 mm): - 11 levels without roller shutter door - 12 levels without roller shutter door H 100 with walls for drawers telescopic: 10 drawers H 80 mm	H 120 with standard dividers (40mm): - 25 levels for trolley with roller - 26 levels for trolley without roller shutter door H 120 with walls for drawers telescopic: 13 drawers H 80 mm
Locking system	Code lock (with spare key) / Badge lock		

Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection
Gross weight (trolley + packaging)	H. 80 : 39 kg H. 100 : 36 kg H. 120 : 38 kg
Dimensions	H. 80 : L. 0,80 x P. 1,20 x H. 1,20 m H. 100 : L. 0,80 x P. 1,20 x H. 1,45 m H. 120 : L. 0,80 x P. 1,20 x H. 1,60 m
Volume	1,15 m ³