

# Neop | 600 x 400 distribution trolleys

## Daily distribution of medicines with pill dispensers and nominative trays

- Robustness: excellent resistance, one-piece structure made of steel (height 80 cm) and of aluminium (height 100 and 120 cm)
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Light
- Ergonomic and handy: designed for medical staff comfort
- Modular: wide range of drawers and accessories
- Customizable: ward name identification



Robust



Hygienic



Ergonomic



Customizable

- Single column trolley
- 3 height: 80, 100 or 120 cm
- Roller shutter door or door wing
- Code lock
- 4 available inside panel models (slides) according to internal fittings

EasyMoove® option: trolley motorization

### Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue



## Standard configurations - PREMIUM

Structure WITH internal fittings and peripheral accessories.  
Add the color reference of the trolley.

### Trolleys height 100 cm

#### DIN – 18 patients bins H. 70

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 18 bins H. 70 l.1/3 – white color
- 7 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



#### Reference

8YM64H100DIN18

#### DIN – 24 patients bins H. 110

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 24 bins H. 110 L. 1/6 – white color
- 5 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 75 mm
- 1 drawer 600 x 400 H. 115 mm
- 1 rail (bottom right)
- 1 bag holder



#### Reference

8YM64H100DI24E1

#### DIN – 24 patients bins H. 70

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 24 bins H. 70 L.1/4 – white color
- 7 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 75
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




#### Reference

8YM64H100DIN24

### DIN – 24 patients bins H. 150

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 1 drawer 600 x 400 H. 75 mm for 60 cups
- 50 transparent cups with lid
- 50 blue cups with lid
- 24 bins L. 1/6 H. 150 mm
- 5 shelves 600 x 400 for bins
- 1 rail (bottom right)
- 1 bag holder



#### Reference

8YM64H100DI24E2

### DIN – 28 patients bins H. 70

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 28 bins H. 70 L. 1/4 – white color
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




#### Reference

8YM64H100DIN28

### PDA PILI – 70 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 batches of 10 PILI dispenser supports
- 7 trays 600 x 400 for PILI delivered without dispenser
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




#### Reference

8YM64H100PIL70

### PDA MODULO CLASSIC – 70 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 7 trays 600 x 400 with 10 Modulo Classic dispenser
- 1 drawer 600 x 400 H115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



#### Reference

8YM64H100MC70

### PDA ROLL DISPENSER – 24 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 24 bins H110 L.1/6 + roll dispenser + divider
- 5 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H75
- 1 drawer 600 x 400 H115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




#### Reference

8YM64H100DEV24

### PDA ROLL DISPENSER – 30 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 30 bins H110 L.1/6 + roll dispenser + divider
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



#### Reference

8YM64H100DEV30

### PDA ROLL DISPENSER – 42 patients

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 42 bins H. 70 L. 1/6 + roll dispenser + divider
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



#### Reference


8YM64H100DEV42

### Trolleys height 120 cm

#### DIN – 36 patients bins H. 110

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 1 drawer 600 x 400 H. 75 mm for 60 cups
- 50 transparent cups with lid
- 50 blue cups with lid
- 36 bins L. 1/6 H. 110 mm – white color
- 7 tablettes 600 x 400 pour bacs tiroirs
- 1 tiroir 600 x 400 H. 155 mm avec séparateurs
- 1 rail (bottom right)
- 1 bag holder




#### Reference

8YM64H120DI36E1

#### PDA ROLL DISPENSER – 36 patients

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 36 bins H. 110 L. 1/6 + roll dispenser + divider
- 7 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket




#### Reference

8YM64H120DEV36

## PDA ROLL DISPENSER – 42 patients

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 42 bins H. 110 L. 1/6 roll dispenser + divider
- 8 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



### Reference

8YM64H120DEV42

## PDA ROLL DISPENSER – 54 patients

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 54 bins H. 70 L. 1/6 roll dispenser + divider
- 10 shelves 600 x 400 for bins
- 1 drawer 600 x 400 H. 115
- 3 rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 2 bag holders
- 1 large side basket



### Reference

8YM64H120DEV54

## Standard configurations – BASIC

Structure WITHOUT internal fittings or peripheral accessories.

Add the color reference of the trolley.

### Trolleys height 80 cm


Reference	Description
8YCHM164H80FR	NEOP Basic medicine trolley 600x400 H80 thermoplastic upper shelf + curtain door + standard dividers (40 mm level)
8YCHM164H80CFR	NEOP Basic medicine trolley 600x400 H80 compact upper shelf + curtain door + standard dividers (40 mm level)
8YCHM164H80FRC	NEOP Basic medicine trolley 600x400 H80 thermoplastic upper shelf + curtain door + code lock + standard dividers (40 mm level)
8YCHM164H80CFRC	NEOP Basic medicine trolley 600x400 H80 compact upper shelf + curtain door + code lock + standard dividers (40 mm level)

## Trolleys height 100 cm

### TROLLEY WITH MANUAL CURTAIN DOOR

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 3 rails (top right, bottom right, bottom left)




Reference	Description
8YCHM164H100FRC	NEOP Basic medicine trolley 600x400 H100 P40 + curtain door + code lock + 3 rails + standard dividers (40 mm level)
8YCHMH100FRC_ED	NEOP Basic medicine trolley 600x400 H100 P62 + curtain door + code lock + 3 rails + space-optimized dividers (62 mm level)

### TROLLEY WITH EASYMOOVE® MOTORIZATION SYSTEM

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 3 rails (middle left, bottom left, bottom right)
- Motorization system




Reference	Description
8YCHM164H100_RCM	NEOP Basic medicine trolley 600x400 H100 P40 + motorization system + curtain door + code lock + 3 rails + standard dividers (40 mm level)
8YCHM164H100P6_RCM	NEOP Basic medicine trolley 600x400 H100 P40 + motorization system + curtain door + code lock + 3 rails + space-optimized dividers (62 mm level)

### TROLLEY WITH ERGONOMIC HANDLE

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 1 ergonomic handle
- 3 rails (middle left, bottom right, bottom left)




Reference	Description
8YCHM164H100_RCP	NEOP Basic medicine trolley 600x400 H100 P40 + curtain door + code lock + 3 rails + ergonomic handle + standard dividers (40 mm level)
8YCHM164H100P6_RCP	NEOP Basic medicine trolley 600x400 H100 P62 + curtain door + code lock + 3 rails + ergonomic handle + space-optimized dividers (62 mm level)

## TROLLEY WITH SLIDING SHELF

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 rails (bottom right, bottom left)
- 1 sliding extension shelf



Reference	Description
8YCHM164H100RCTC	NEOP Basic medicine trolley 600x400 H100 P40 + curtain door + code lock + sliding shelf + 2 rails + standard dividers (40 mm level)
8YCMH100FRCTC_ED	NEOP Basic medicine trolley 600x400 H100 P62 + curtain door + code lock + sliding shelf + 2 rails + space-optimized dividers (62 mm level)

## TROLLEY WITH FOLDABLE SHELF

600 x 400 medicine trolley H. 100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 rails (bottom right, bottom left)
- 1 foldable extension shelf



Reference	Description
8YCHM164H100RCT	NEOP Basic medicine trolley 600x400 H100 P40 + curtain door + code lock + foldable shelf + 2 rails + standard dividers (40 mm level)
8YCMH100FRCT_ED	NEOP Basic medicine trolley 600x400 H100 P62 + curtain door + code lock + + foldable shelf + 2 rails + space-optimized dividers (62 mm level)


## Trolleys height 120 cm

### TROLLEY WITH COMPACT LAMINATE TOP SHELF

600 x 400 medicine trolley H. 120 equipped with:

- Compact laminate top shelf
- Roller shutter door
- Code lock

Reference	Description
8Y164H120P6C_RC	NEOP Basic medicine trolley 600x400 H120 P62 + curtain door + code lock + compact top shelf



### TROLLEY WITH THERMOPLASTIC UPPER SHELF

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 3 rails (top right, bottom right, bottom left)

Reference	Description
8YCHM164H120FRC	NEOP Basic medicine trolley 600x400 H120 P40 + curtain door + code lock + 3 rails
8YCHMH120FRC_ED	NEOP Basic medicine trolley 600x400 H120 P62 + curtain door + code lock + 3 rails




### TROLLEY WITH FOLDABLE SHELF

600 x 400 medicine trolley H. 120 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 2 rails (bottom right, bottom left)
- 1 foldable extension shelf

Reference	Description
8YCHM164H120RCT	NEOP Basic medicine trolley 600x400 H120 P40 + curtain door + code lock + foldable shelf + 2 rails
8YCMH120FRCT_ED	NEOP Basic medicine trolley 600x400 H120 P62 + curtain door + code lock + foldable shelf + 2 rails






# Custom configurations

## STEP 1 - Trolley structure

This reference includes structure, wheels and inside panels.

	Trolley height 80 cm*	Trolley height 100 cm*	Trolley height 120 cm*
Standard inside panels for standard drawers (40 mm level)	8Y164H80P4	8Y164H100P4	8Y164H120P4
Inside panels for telescopic drawers	-	8Y164H100T	8Y164H120T
Inside panels for « Adhesia Equipments » pill dispensers (40 mm level)	-	8Y164H100W	8Y164H120W

## STEP 2 - Trolley customization

		Trolley height 80 cm*	Trolley height 100 cm*	Trolley height 120 cm*	
1	Upper shelf	Thermoplastic shelf with cutout	8Y001_6400TH		
		Thermoplastic shelf with cutout compatible with sliding shelf	-	8Y001_6400TC	-
		Thermoplastic shelf without cutout	8Y001_6400THSP		
		Thermoplastic shelf without cutout compatible with sliding shelf	-	8Y001_6400TCSP	-
		Compact laminate shelf with cutout	8Y001_6400CP1		
2	Closing system	Manual curtain door	8Y002_6008F3	8Y002_6010F3	8Y002_6012F3
		Motorized curtain door Neop® Lock	-	9U000	-
		Door wing for telescopic drawers	-	8Y002_0010F2	8Y002_0010F2
		Door wing for standard dividers	8Y002_0008F1	8Y002_0010F1	-
3	Code lock	For curtain door	8Y003CODE3		
		For door wing and telescopic drawers	8Y003CODE2		
		For door wing and standard dividers	8Y003CODE1		
4	Colour	Magenta pink 	8Y004_0008RM	8Y004_0010RM	8Y004_0012RM
		Apple green 	8Y004_0008VP	8Y004_0010VP	8Y004_0012VP
		Sky blue 	8Y004_0008BC	8Y004_0010BC	8Y004_0012BC
5	Specific marking	8Y132			
6	EasyMoove® trolley motorization	-	9R002_L01	9R002_L01	

\*For precise height of working surface, see technical data


## Options, accessories and internal fittings

### EasyMoove® trolley motorization

- 100% motorized driving
- Limits pushing efforts for users
- Reduces the risk of musculoskeletal disorders



### Extension shelf, ergonomic handle, IT support, upper frame, drawers...



# Technical data

## Trolley with thermoformed plastic top shelf

### Overall dimensions of the trolley



### Useful dimensions of the worktop



*Ergonomic gripping area*

### Worktop height

H. 120 : 1297 mm  
H. 100 : 1057 mm  
H. 80 : 897 mm

## Trolley with compact laminate top shelf

### Overall dimensions of the trolley



### Useful dimensions of the worktop



### Worktop height

H. 120 : 1304 mm  
H. 100 : 1064 mm  
H. 80 : 904 mm

The overall dimensions are taken with the wheels in the least bulky position.  
The depth includes a lock by default.

## Materials

<b>Structure</b>	20/10th complete welded one-piece structure, Epoxy anti-corrosion treated H 80: Steel (electro-galvanized) / H 100: Aluminium / H 120: Aluminium
<b>Work surface</b>	Thermoplastic: ABS / PMMA OR compact laminate th. 10 mm
<b>Slides</b>	ABS (Acrylonitrile Butadiene Styrene) OR ballbearing slide made of galvanized steel with shock absorber
<b>Closing system</b>	Roller shutter door: PP (Polypropylene), 2 mm thickness Door wing: Aluminium with EPOXY painting
<b>Wheels</b>	Simple non marking rubber wheels - Polyamide cover (PA)
<b>Bumpers</b>	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)
<b>Strip of colour</b>	70 micron calendered PVC vinyl film, with acrylic adhesive

## Additional information

<b>Wheels</b>	125 mm diameter - 2 breaking wheels + 2 swivelling wheels - wheel and wheel-pin double blocking		
<b>Empty weight</b>	H. 80: 27 kg	H. 100: 26,5 kg	H. 120: 28,8 kg
<b>Capacity</b>	<b>H 80 with standard dividers (40 mm):</b> - 15 levels of 40 mm for trolley with roller shutter door - 16 levels of 40 mm for trolley without roller shutter door	<b>H 100 with standard dividers (40 mm):</b> - 19 levels for trolley with roller shutter - 20 levels without roller shutter door  <b>H 100 with dividers at optimized size (62 mm):</b> - 11 levels without roller shutter door - 12 levels without roller shutter door  <b>H 100 with walls for drawers telescopic:</b> 10 drawers H 80 mm	<b>H 120 with standard dividers (40 mm):</b> - 25 levels for trolley with roller - 26 levels for trolley without roller shutter door  <b>H 120 with walls for drawers telescopic:</b> 13 drawers H 80 mm
<b>Locking system</b>	Code lock (with spare key)		

## Packaging

<b>Packing</b>	Individually on pallet strapped and film-wrapped, OR on wheels, with proper protection
<b>Gross weight (trolley + packaging)</b>	H. 80: 39 kg H. 100: 38,5 kg H. 120: 40,8 kg
<b>Dimensions</b>	H. 80: L. 0,80 x P. 1,20 x H. 1,20 m H. 100: L. 0,80 x P. 1,20 x H. 1,45 m H. 120: L. 0,80 x P. 1,20 x H. 1,60 m
<b>Volume</b>	1,15 m <sup>3</sup>