

Clean linen range

Nursing trolleys

Photos are not contractual

The range of nursing trolleys streamlines and facilitates the work of nursing assistants. The **ergonomic** design and **modular** concept of storage compartments make it easy to adapt the trolley to each ward's needs.

To limit the pushing effort required from users and prevent musculoskeletal disorders, our trolleys can be equipped with the **EasyMoove motorization system**. **Computerized management** is possible by adding computer or tablet holder to the trolley.



Advantages

- ▶ Open or closed model with sliding doors
- ▶ 2 models: 15 to 20 patients and 20 to 30 patients
- ▶ Adjustable-height middle tray with removable dividers
- ▶ Several storage areas
- ▶ Materials resistant to corrosion and disinfecting products



Hygienic



Ergonomic

Clean linen range

Nursing trolleys

Mini nursing trolleys: 15 to 20 patients

Standard equipments

- ▶ 1 300 x 400 height 100 sliding ISO tray under upper shelf
- ▶ 1 adjustable and removable middle shelf
- ▶ 2 removable dividers for the middle shelf
- ▶ 4 standard rails
- ▶ 1 storage bin
- ▶ 2 large side baskets

Description	Reference
Mini nursing trolley without doors	8P500
Mini nursing trolley without doors with upper divider	8P500+
Mini nursing trolley without doors with EasyMoove® motorization system	8P500_EM
Mini nursing trolley with doors and code lock	8P500_C
Mini nursing trolley with doors, code lock and upper divider	8P500+_C
Mini nursing trolley without doors with code lock and EasyMoove® motorization system	8P500_CEM
Mini nursing trolley with doors and key lock	8P500_P
Mini nursing trolley with doors, key lock and upper divider	8P500+_P
Mini nursing trolley with doors, key lock and EasyMoove® motorization system	8P500_PEM

Dimensions



Clean linen range

Nursing trolleys

Classic nursing trolleys: 20 to 30 patients

Standard equipments

- ▶ 1 300 x 400 height 100 sliding ISO tray under upper shelf
- ▶ 1 adjustable and removable middle shelf
- ▶ 3 removable dividers for the middle shelf
- ▶ 4 standard rails
- ▶ 1 storage bin
- ▶ 2 large side baskets

Description

Reference

Classic nursing trolley without doors	8P550
Classic nursing trolley without doors with upper divider	8P550+
Classic nursing trolley without doors with EasyMoove® motorization system	8P550_EM
Classic nursing trolley with doors and code lock	8P550_C
Classic nursing trolley with doors, code lock and upper divider	8P550+_C
Classic nursing trolley with doors, code lock and EasyMoove® motorization system	8P550_CEM
Classic nursing trolley with doors and key lock	8P550_P
Classic nursing trolley with doors, key lock and upper divider	8P550+_P
Classic nursing trolley with doors, key lock and EasyMoove® motorization system	8P550_PEM

Dimensions



Clean linen range

Nursing trolleys

Technical data

Materials

Structure	Structure: 40x20 steel tube - Shelves: 1.5mm thick welded steel sheet The entire assembly is treated with epoxy anti-corrosion coating
Wheels	Dual-band Polyamide (PA) wheels - Flexible Polyurethane (PUR) tread
Bumpers	4 round bumpers made of non-marking EVA (Ethylene Vinyl Acetate)
Doors	Translucent PMMA, 6 mm thick
Removable shelf	Steel wire Ø 6 and Ø 4 treated with epoxy anti-corrosion coating
Adjustable dividers	6 mm diameter steel wire with epoxy anti-corrosion treatment
Storage bin	ABS (Acrylonitrile Butadiene Styrene)
Drawer	ABS injection-molded drawer 300 x 400
Basket	Steel coated with epoxy polyester RAL7024

Additional information

Wheels	4 swivelling wheels Ø 125 mm (1 with brake and 1 with directional lock)			
Pre-equipment	To be struck on the plateau			
Empty weight	8P500: 38 kg 8P550: 51 kg	8P500+: 41 kg 8P550+: 57 kg	8P500_P: 51 kg 8P550_P: 65 kg	8P500+_P: 54 kg 8P550+_P: 68 kg
Load capacity	Top shelf: 60 kg - Removable shelf: 30 kg - Lower shelf: 60 kg			
Locking system	For trolleys with doors: key lock or code lock (with spare key)			

Packaging

Packing	Individually on pallet + cardboard packaging	
Gross weight	8P500: 50 kg	8P550: 66 kg
Dimensions	L. 0,80 x D. 1,20 x H. 1,30 m	
Volume	1,3 m³	

Colors

Description		Reference
Pop blue color		8P504_BP
Magenta pink color		8P504_RM
Blueberry color		8P504_MY
Orange color		8P504_OR
Apple green color		8P504_VP

Clean linen range

Nursing trolleys

Motorization


All our nursing trolleys can be equipped with the motorized wheel EasyMoove®. This innovation enables to **reduce physical efforts of the users**. A simple act with the finger on the handle allows the system to operate. Staff's work is therefore with utmost convenience.

Description	Reference
EasyMoove® motorization system for Mini nursing trolley (to be ordered at the same time as the trolley)	9R002_L02V2
EasyMoove® motorization system for Classic nursing trolley (to be ordered at the same time as the trolley)	9R002_L03V2



IT supports

Description	Reference
1 Sheet metal laptop holder	8Y009CP4
2 PMMA laptop holder with cover	8Y009CP2
3 PMMA laptop holder without cover	8Y009CP1
4 Touch tablet holder	8Y009CP3
5 Vesa screen holder	8Y009CP5
6 Computer arm extension	QUI2F110
Compact mouse holder for compact laptop holder	28502
Cable guide	8E105
10 cable ties	8E109





Clean linen range

Nursing trolleys

Accessories for standard rails

	Description	Overall dimensions (in mm)	Reference
	Side basket small model	L. 355 x P. 218 x H. 238	28062N.GR
7	Side basket big model	L. 405 x P. 218 x H. 238	28061N.GR
8	Storage bin	L. 455 x P. 135 x H. 110	28063N.GR
9	Chrome wire briefcase		BF20250
10	Steel document holder	L. 352 x P. 91 x H. 160	28013
11	Vertical support for allocation sheets		AN09403
12	Waste bin holder	L. 173 x P. 265 x H. 302	28080.GR
13	Pedal waste bin 25 L	L. 235 x P. 190 x H. 400	28074N
14	Standard bracket for gloves	L. 135 x P. 80 x H. 300	28001N.GR
15	Removable roll support		28059N.GR



Autres équipements

	Désignation	Référence
16	Handle	8P502
17	Removable 3-compartment upper divider	8P515
18	Bottle holder - small model	28038
	Removable shelf for Mini nursing trolley	ACC8P505
	Removable shelf for Classic nursing trolley	ACC8P554
	Fixed bag holder with manually opening lid	8P525
19	Porte-sac fixe avec couvercle à ouverture manuelle	8P525
20	Porte-sac amovible avec couvercle et pédale d'ouverture	22001_S
	White fabric laundry bag	22100_S

