

Compression range

Pinel® restraint system

The **Pinel®** design allows you to address any known restraint requirement with ONE system. At the emergency end of restraining, the aggressive patient can be secured to a seven-point restraint level in about 10 seconds ; but, adjusted safely to a comfortable one-point restraint as he/she calms. At the passive end of restraint needs, confused patients can be secured with full mobility within the bed or chair to the level that patients are not even aware of being restrained.

This single system addresses all possible restraint needs within or between (transport) any institution, e.g. Long Term Care, Emergency Room, Psychiatric, ICU or Correctional facilities...



Advantages

- ▶ Confort: Polyester + Kevlar
- ▶ Security: magnetic locking system
- ▶ Hygiene: treatment against bacteria, easy washing & fast drying
- ▶ Liberty of movements: very flexible when the patient is secured
- ▶ Rapidity: universal system

Compression range

Pinel® restraint system

Complet kit

The kit contains:

- ▶ 4 limb straps
- ▶ 1 waist belt
- ▶ 2 side straps
- ▶ 1 extender
- ▶ 1 belt for trunk
- ▶ 1 pelvic strap
- ▶ 1 strap belt
- ▶ 3 keys
- ▶ 13 axis + bolt with string
- ▶ 2 axis + bolt
- ▶ 1 emergency cutter

Equipped with a color code, this case permits you to put easily the different equipments of the Pinel® range.

Description	Reference
Complet kit with case	SAA50100
Complet kit without case	SAA50200
Case only	SAA50900

Bolts and accessories

Every Pinel® fixation is equipped with very effective locking system which can quickly be applied and easily used. The locking system with axis + bolt enables to block until 4 straps together. The system unlocks by putting the magnetic key in the bolt.

Resistance: 500 kg.

Description	Reference
Magnetic key	SAA50910
Axis + bolt	SAA50940
Axis + bolt + key	SAA50950

Pelvic strap



This abdominal belt allows the patient to lateralize and move up and down. For patients who are in danger of sliding in a bed or on a chair.

The kit contains 1 abdominal belt, 2 free side straps, 1 pelvic strap, 2 axis + bolt with string, 2 axis + bolt and 3 keys.

Description	Reference
Abdominal belt kit small size (from 38 to 97 cm)	SAA50402
Abdominal belt kit standard size (from 53 to 112 cm)	SAA50400
Abdominal belt kit large size (from 78 to 137 cm)	SAA50401
Extension with abdominal belt (recommended for the obese - L. 30 cm)	SAA50410
Pelvic strap	SAA50420

Limb strap



The kit contains 2 limb strap, 1 magnetic key et 4 axis + magnetic bolt and twine. The limb strap are equipped with a velcro closure system allowing quick and easy control wrists and ankles in an emergency situation.

Description	Reference
Limb strap kit small size	SAA50301
Limb strap kit standard version	SAA50300