

Patient file management

XXL medical file trolley

The XXL patient file trolley range is designed for the storage of medical and x-ray files. These ergonomic and modular items enable an optimised storage of the files according to each ward's needs. File confidentiality is assured thanks to the curtain door closing.

You can also adapt a PC support on all trolleys and therefore get access to computerized patient files.

Non contractual photos



Hygienic



Ergonomic

Advantages

- ▶ Adapted to medical files with large size
- ▶ Ergonomic and modular
- ▶ Guaranteed confidentiality
- ▶ Patient file computerisation
- ▶ Resistant to the disinfecting products of hospital environment

XXL medical file trolley


Trolley

Description

XXL medical file trolley

Reference

8C010BP



Technical data

Materials

Structure	Anodized aluminum profile
Side panels	White PS (Polystyrene)
Work surface	ABS (Acrylonitrile Butadiene Styrene)
Curtain door	Aluminum end blade - PP (Polypropylene) plastic blade
Handle	Aluminum tube Ø 20 - Nylon end cap
Wheels	Simple bandage in a Polyamide (PA) cap
Bumpers	4 non-marking EVA (Ethylene Vinyl Acetate) bumpers
Partitions	Recycled PS (Polystyrene)

Additional information

Wheels	4 wheels Ø 125 mm - 2 with brakes and 2 swiveling
Locking system	Key lock
Empty weight	61 kg
Maximum capacity	2 x 33 files of 20 mm 2 x 16 files of 40 mm 2 x 11 files of 60 mm (standard) 2 x 8 files of 80 mm 2 x 7 files of 100 mm
Partitions	Adjustable partitions in 20 mm increments Number of increments per level: 33 Number of partitions supplied: 20

Packaging

Packing	Individually on pallet strapped and film-wrapped
Gross weight (trolley + packaging)	71 kg
Dimensions	L. 0,80 x D. 1,20 x H. 1,40 m
Volume	1,3 m ³




XXL medical file trolley

Optional equipment

Rail and compatible accessories

These options can only be installed on the right side of the trolley.


	Description	Overall dimensions (in mm)	Reference
1	Standard rail		8C009
2	Small wire basket	L. 355 x D. 218 x H. 238	28062N.GR
	Big wire basket	L. 405 x D. 218 x H. 238	28061N.GR
3	Storage bin	L. 455 x D. 135 x H. 110	28063N.GR
4	Document holder	L. 350 x D. 106 x H. 160	28013



IT supports

These options can only be installed on the right side of the trolley.

	Description	Reference
5	Sheet metal laptop holder	2F100
6	PMMA laptop holder with cover	2F105
7	PMMA laptop holder without cover	2F101
8	Touchscreen tablet holder	2F103
9	VESA screen holder	2F106
10	Compact laminate mouse pad for compact laptop holder	28502



Other equipment

Description	Reference
Removable partition	ACC8C020

