

Accessories and equipment

Neop® sliding shelf

Non-contractual photos

The sliding shelf is an essential element on a care and medication trolley. It offers an additional work surface, adapting to the changing needs of a medical environment.

The shelf can be deployed to the right or left, thus optimizing the use of the trolley's equipment. When not in use, it easily retracts to free up space.

The Neop® sliding shelf contributes significantly to the ergonomics of the caregivers' workstation, thus facilitating their daily tasks.



Hygienic



Ergonomic

Benefits

- ▶ Fully integrated under the worktop
- ▶ Opens to the right or left
- ▶ Easy to deploy and fold thanks to the castor guidance mechanism
- ▶ Robust: designed to support loads suitable for hospital environments
- ▶ Easy to clean surface without roughness, resistant to scratches and daily wear
- ▶ Does not compromise the functionality or stability of the trolley

Neop® sliding shelf

Technical features

- ▶ Dimensions:
 - Depth: 390 mm
 - Output width: 407 mm
 - Thickness: 10 mm
- ▶ Main material: white smooth compact tablet
- ▶ Assembly type: slide on guide roller for easy deployment and folding
- ▶ Deployment mechanism:
 - Sliding system left or right with automatic locking
 - Manual opening / closing
- ▶ Load capacity: 5 kg
- ▶ Fixing: on the underside of the trolley worktop using a stop system
- ▶ Compatibility: Neop® trolleys 600 x 400 H. 100 cm (this option raises the trolley by 15 mm)

Description

Reference

Neop® sliding shelf

8Y008_0410

